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The irradiated surface will become warm after continuous long-term irradiation. Rest for about **1minute** until it is sufficiently cooled.

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Irradiation to the same area & site may cause redness due to excessive heat buildup on the skin. Let's remove hair with **one round trip**.

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Make sure that the irradiated surface is in **close contact with the skin**. Floating from the skin will be easier to feel pain and heat.

④ 入れ墨・タトゥー等への照射厳禁 Irradiation of tattoos is strictly prohibited

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Never irradiate tattoos as they may cause burns. Please note that we cannot take responsibility.

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To customers after treatment

The Precautions After Hair Removal

- The skin will be sensitive after hair removal, pores are filled with redness and heat.
- Warm up body to improve the blood circulation, you may feel the redness and an itch easily.
- Please keep "aftercare" in your mind even if you were already at home.

Pay attention to the dryness, irritation and friction.

AFTERCARE

Cautions for 3 days since the operation day

- Keep your skin moisturized so as to make the 'hair removal effect UP'.
- Avoid taking shower on the day of hair removal try to take warm bath and do not rub your body too hard.
- Avoid sauna, bedrock bath and sports which may make you sweat.
- Pay attention that the sunburned skin will be easy to receive the stimulus.
- Keep away from drinking alcohols and using the razor.

When you feel the redness, an itch or the tingling, please use the ice or the cold water to cool down.
If there are any other matters weigh on your mind, please contact the salon.

The pace of hair removal

1 time every 2 weeks for half a year.

There are individual differences in terms of hair removal effect.

However, the effect could be expected if the hair removal machine is used in right way.

First please do 'Self Hair Removal' 1 time every 2 weeks for half a year.

Then putting concentration on the parts you are concerned about would be recommended.

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